Kit for semi-quantitative detection of contaminant microorganisms in must, beer or water

Developed and patented by BiorSelf srl

Description:


Self-Beer® is an innovative device for semi-quantitative detection of contaminant microorganisms in must and beer. It can determine the presence of wild yeasts, mould, lactic acid and acetic acid bacteria. Self-Beer® is a union between a species-specific culture-dependent detection method and simple ready-to-use technology (dip-slide) that allows the producer to perform a **rapid self-controlled analysis**.

Advantages:

1. **Autonomy:** Self-Beer® tests can be carried out directly on the production site without requiring laboratories or specialized individuals, while obtaining a clear, evidence-based result
2. **Clear results:** the user can evaluate the presence of contaminant microorganisms in beer through observing colonies grown in the medium
3. **Price:** low costs compared to laboratory analysis

Composition of the kit and uses:


Evaluation of the presence of contamination by wild yeasts, lactic acid bacteria and acetic bacteria. It is also possible to assess the viability of *Saccharomyces* yeasts. The self-beer® analysis kit consists of two test tubes, one green and one red. The tube with the green cap allows the analysis of yeasts: on side A, total yeasts grow, on side B, non-*Saccharomyces* yeasts (wild yeasts) grow. The test tube with the red cap allows the analysis of bacteria: on side A, acetic bacteria grow, while on side B, lactic acid bacteria grow.

Instructions and recommendations for use:

- unscrew and extract the slides, avoiding contact with the agar surfaces;
- Fill the containers with the liquid to be analysed until the line on the label;
- Shake gently and wait 60 seconds;
- Empty and drain Self-Beer® on a piece of paper;
- Close and leave at room temperature (18-30 °C) for 3 – 5 days;
- observe the result obtained (no need for microscopes or magnifying devices)

Instructions for use

Self-beer® is designed to be used in the brewery under normal operating conditions. However, a series of precautions and good sampling practices are recommended:

- sanitize the sampling nozzles;
- purge a little of product before filling Self-beer®;
- if necessary, in some stages of sampling, use well sanitized jugs or containers;
- avoid keeping Self-beer® open in the air for longer times than those necessary for its filling and emptying;
- do not touch the culture media or the inside of the tubes with hands

A non optimal choice of sampling phase can lead to dubious results, in particular, sampling during the phase of full fermentation, whether in tank or bottle, should be avoided. In the case of refermentation, sampling should be carried out immediately before the addition of yeast or when the product is stabilised. For more information, please contact us.

Tel +39 3518614106
E-mail info@biorself.it

BiorSelf Srl
Analysis by Yourself
P.Iva: 03784790788


www.biorself.it

Results interpretation:

If the tested sample is contaminated, after 3-5 days at room temperature (18-30 °C), it will be possible to observe colonies (visible as dots of various sizes and colours) on the surface of the agarised medium. On side A of the green test tube and side B of the red test tube, it will also be possible to see a colour change to yellow, in case of microbial growth.

Evaluate contamination:

If the sample is contaminated, the number of colonies present in the medium will indicate a semi-quantitative measure of the contamination level.

LOW CONTAMINATION	MEDIUM CONTAMINATION	HIGH CONTAMINATION
LOW	10 – 10 ² CFU/ml	
MEDIUM	10 ³ – 10 ⁴ CFU/ml	
HIGH	≥ 10 ⁵ CFU/ml	

Storage:

Store in 10-25°C conditions, away from light, until the expiration date on the external side of box.

Shelf-life:

The 'best before date' is set at **12 months** from the date of production. As the product is an agarised gel, it may be subject to evaporation and natural weight loss depending on seasonal weather and storage/transport conditions.

Quality control:

Each Self-Beer® batch was tested with the following microorganism strains:

Strain Cod	Species	Growth
DSM 70451	<i>Saccharomyces cerevisiae</i>	Green lid side A
DSM 2548	<i>Saccharomyces pastorianus</i>	Green lid side A
DSM 21959	<i>Pichia membranifaciens</i>	Green lid side B
DSM 2002	<i>Acetobacter aceti</i>	Red lid side A
DSM 20019	<i>Lactobacillus curvatus</i>	Red lid side B

Frequently asked questions (F.A.Q.):

- The colour of the culture media has changed, is this normal? Yes, the colour of the culture media for total yeast (green cap side A) and lactic acid bacteria (red cap side B) will change as the microorganisms grow.
- Is it possible to assess the viability of the yeast inoculum for fermentation with self-beer®? Yes, *Saccharomyces cerevisiae* yeasts grow in the medium on side A of the green lid tube. If the yeast is viable and uncontaminated, you should get a very high growth on side A and no colony on side B (wild yeast medium).
- The colour has changed but I cannot see the colonies, why is that? In the case of very high contamination (such as immediately after yeast inoculation) it is possible that the colonies may form a uniform layer and thus make the individual colonies not visible. In this case, an opaque patina can be observed.

Warnings and precautions:

The product must be kept out of the reach of children and must be used according to the manufacturer's documentation and specifications.

Self-beer® is not classifiable as hazardous, please read the documents inside the packaging for more information.

Do not use Self-Beer® with samples at high temperatures (>50 °C) to avoid deterioration of the media.

BiorSelf will not be responsible in any way for uses that do not conform to these specifications. It is the responsibility of customers to act in accordance with applicable health, environmental and safety regulations and to take appropriate action in relation to the storage, handling, sale and use of products.

Disposal of waste:

Disposal of the product must be carried out in accordance with current national and local regulations.

The product is not classifiable as special hazardous waste according to Regulation (EU) 1357/2014.

The case must be disposed of in the paper collection - PAP21.

The tubes must be disposed of in the unsorted municipal collection.

Tel +39 3518614106
E-mail info@biorself.it

BiorSelf Srl
Analysis by Yourself
P.Iva: 03784790788

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Tel +39 3518614106
E-mail info@biorself.it

BiorSelf Srl
Analysis by Yourself
Piva: 03784790788

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